

# SMARTACT SMARTACT PIN




The introduction of the Smartact mechanical handpiece and the Smartact Pin titanium tacks revolutionise the procedures and the traditional manual instruments used in regenerative oral surgery.

## Smartact

The unique characteristics of Smartact enable the clinician to improve and simplify the aspects regarding fixing the membrane in the precision of the surgical operation, in access to all areas of the oral cavity and in the reduction of surgical times.

## Smartact Pin

The Smartact Pin design and connection enable axial insertion of the tack due to the single-use stop located on the handpiece that facilitates fixing the membrane in any situation and clinical area, offering the best possible stability and protection of the graft.



courtesy of Dr. Ferdinando D'Avenia



## Advantages of the Smartact instrument

- Mechanical handpiece
- Pedal action
- Ergonomic and autoclavable tip
- Excellent access to all the oral cavity areas
- Simple to Use and surgical times reduced
- Maximum precision in determining fixing point of the membrane
- Adjustment of the positioning force of the tack in the various bone densities
- Greater patient comfort
- Reduction of risk factors

**Ref. 4850** Smartact handpiece and pedal  
**Ref. 4730** Smartact Pin 3 titanium tacks L 3 mm and  
1 single-patient positioner - Sterile.



Art and Innovation in Medical Technology

CGM spa  
**DIVISIONE MEDICALE META**  
VIA E. VILLA 7 42124 REGGIO EMILIA ITALY  
TEL. +39 0522 502311 FAX. +39 0522 502333  
www.metahosp.com info@metahosp.com

# SMARTACT SMARTACT PIN



**Innovation,  
Simplification and Safety  
in the fixing  
of membranes.**

**Improved surgical  
performance  
and absolute comfort  
for the patient.**

